




## Compact Dry™ LS

Ready-to-Use Medium for  
*Listeria ssp.*



### Background

It is important to detect *Listeria ssp.* in food products and the food environment to limit the possibility of food poisoning. A mixing and dilution culture method has been widely used to determine microbial count. The method is time-consuming and requires complicated operations such as preparation of hot agar, mixing a dilution uniformly and/or spreading. To save operator time and make it possible for anyone to perform a bacterial count test without difficulty, Compact Dry was developed based on a new concept and technology applicable to the food industry.

### Detection

Compact Dry LS detects *Listeria ssp.* in both food and environmental samples. Results are measurable after a 24-hour incubation.

*Listeria ssp.* appear as light blue/blue colonies on the plate.

Colonies sampled from the Compact Dry plate can be used for confirmation of *Listeria ssp.* by inoculation of colonies onto selective media.

### Warnings and Precautions

#### 1. General precautions

- Read and follow precisely the warnings and directions for use described in the package insert and/or label.
- Do not use the product after its expiration date. Quality of the product is not warranted after its shelf life expires.
- Do not use product that contains any foreign materials, is discolored or dehydrated, or has a damaged container.
- Use plates as soon as possible after opening. Return any unused plates to the aluminum bag and seal with tape to avoid light and moisture.
- Cap tightly after inoculation to avoid dehydration of gelled medium.

#### 2. Safety precautions

- If medium or reagent comes into contact with eyes or mouth, immediately wash with water and consult a physician.
- Procedures with microorganisms involve certain risks of laboratory-acquired infections. Procedures should be carried out under the supervision of trained laboratory personnel with biohazard protection measures.
- Treat any laboratory equipment or medium that comes into contact with the specimen as infectious and sterilize appropriately.

#### 3. Precautions for disposal of waste

- Sterilize any medium, reagent or materials by autoclaving or boiling after use, and then dispose of it as industrial waste according to local laws and regulations for disposal of such material.

#### 4. User responsibilities

- It is the user's responsibility in selecting any test method to evaluate a sufficient number of samples with particular foods and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.
- It is the user's responsibility to determine that any test methods and results meet its customers or suppliers' requirements. The user must train its personnel in proper testing techniques.
- It is the user's responsibility to validate the performance of this method for use with any non-certified matrix.

#### 5. Limitation of warranties

- Compact Dry plates are manufactured at ISO 9001:2015 facility. If any Compact Dry plate is proven to be defective by fault of the manufacturer or its authorized distributors, they may replace or, at their discretion, refund the purchase price of any plate. These are the exclusive remedies.

### Storage and Shelf Life

Storage: Keep at room temperature (1–30°C)

Shelf life: Eighteen (18) months after manufacturing. Expiration date is printed on outer box label and aluminum bag label.

### Package

Compact Dry LS 100 plates      Code 54060  
Compact Dry LS 1400 plates      Code 54060-CS

### Further Information

#### Customer Support

Shimadzu Diagnostics Corporation  
3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan  
Phone: +81-3-5846-5707  
contact@sd.c.shimadzu.co.jp

### Manufactured by

Shimadzu Diagnostics Corporation  
3-24-6, Ueno, Taito-ku, Tokyo 110-0005, Japan

Kit components, operating  
instructions and interpretation



## Operating Procedure

### Preparation of specimen

- Bacterial number in solid foods:** Add buffering solution to the sample and homogenize. Drop 1 ml of specimen (to be further diluted if necessary) in the middle of a Compact Dry plate.
- Bacterial number in water or liquids:** Drop 1 ml of specimen (to be diluted if necessary) in the middle of a Compact Dry plate.
- Bacterial number in swab test specimen:** Drop 1 ml of wiping solution (to be diluted if necessary) obtained from a cotton swab, in the middle of a Compact Dry plate. It is recommended to use Easy Wiping Kit, available as an optional kit.
- For effective recovery of *Listeria* spp., it is recommended that samples be inoculated to Compact Dry after 1 hour of resuscitation step at 20°C by Buffered Peptone Water. Solid food samples should be homogenized by BPW. Liquid samples and swab samples should be inoculated to a small amount of BPW. This resuscitation step does not affect the bacterial number according to ISO 11290-2 method for the detection and enumeration of *Listeria monocytogenes*.

**pH Adjustment:** The pH of the product or 1:10 dilution of product should be between 6 and 7 for optimal growth of target microorganisms. If the pH is not between 6 and 7, adjust the pH of the product or 1:10 dilution with 1 N or 0.1 N NaOH for acidic products or 1 N or 0.1 N HCl for alkaline products.

## Directions for Compact Dry LS

- Open aluminum bag, and take out a set of 4 plates.
- Detach necessary number of plate(s) from a set by bending up and down while pressing the lid. Use a set of four connected plates when serial dilution measuring is intended.
- Remove cap, and pipette 1 ml of sample in the middle of dry plate. Replace cap. The specimen diffuses automatically and evenly on the sheet (total medium of 20 cm<sup>2</sup>) to transform it into gel within seconds.
- Turn over the capped plate and incubate 24 ± 2 hours or 48 ± 3 hours at 35 or 37 ± 1°C.
- From the backside of the plate, count the number of light blue/blue colonies in the medium. White paper placed under the plate can make colony counting easier. For large numbers of colonies, use the grids carved on the backside consisting of 1 cm x 1 cm, or 0.5 cm x 0.5 cm, at the four corners.

## Precautions for Use

- During inoculation, do not touch the surface of plate, and be careful to avoid any contamination by falling microorganism.
- During incubation, keep the lid tight to avoid any possible dehydration.
- It is recommended to use a stomacher bag with filter to eliminate risks of contamination with tiny pieces of food.
- Specimen should be diluted by buffer solution to the level of concentration of less than 300 cfu/plate.
- If bacteria more than 10<sup>4</sup> cfu were inoculated on a plate, no colonies are formed, but the entire plate will become colored.
- If the nature of sample affects the reaction of the medium, inoculate the sample only after the factor is eliminated by means such as dilution. For instance, samples with high viscosity, colored, reacted with redox indicator, and too high or too low a pH.

## Interpretation

*Listeria* spp. forms light blue/blue colonies of 1–2 mm in diameter by chromogens contained in the medium.

## Precaution for interpretation

- The full plate size is 20 cm<sup>2</sup>. The plate's backside contains a carved grid of 1 cm x 1 cm and 0.5 cm x 0.5 cm to make colony counting easier. If large numbers of colonies are present on the medium, the total viable count can be obtained by averaging the number of colonies per large grid (1 cm x 1 cm), counted from several grids, and multiplying by 20.
- It is known that *Listeria ivanovii* tends to grow slowly and *Listeria seeligeri* tends to be inhibited in this plate.