

## Centrifuge Tubes

Sterile, highly durable, and suited for a wide temperature range



### Applications

- Cell centrifugation
- Molecular separation
- Light-sensitive sample storage
- Handling of specialized reagents

### Technical Specifications

|                               |                                                                                  |
|-------------------------------|----------------------------------------------------------------------------------|
| <b>Material Type</b>          | PP (tube) + HDPE (cap)                                                           |
| <b>RNase/DNase</b>            | RNase/DNase free                                                                 |
| <b>Pyrogens</b>               | Pyrogen free                                                                     |
| <b>Temperature Range (°C)</b> | -81°C to 121°C                                                                   |
| <b>Surface Properties</b>     | Ultra hydrophobic surface                                                        |
| <b>Centrifugal Strength</b>   | Withstands up to 12,000 × g                                                      |
| <b>Sterilization</b>          | E-beam or ionizing radiation                                                     |
| <b>Printed Graduations</b>    | Easy-to-read printed graduations, accurate to 33%, with dedicated labeling space |
| <b>Shelf Life</b>             | Three-year shelf life from manufacturing date                                    |
| <b>Rack</b>                   | Available with or without rack, made from virgin polypropylene                   |


## Centrifuge Tubes

### Quality Assurance Certification


Products are manufactured in compliance with international ISO 9001:2015, ISO9001, and ISO13485 standards. Free from R-NASE, D-NASE, and pyrogens, meeting USP Class VI requirements.

### Technical Drawing

15mL Centrifuge Tubes / Material meets USP Class VI standards



50mL Centrifuge Tubes / Material meets USP Class VI standards



### Ordering Information

| Code    | Description                    | Packaging Details          |
|---------|--------------------------------|----------------------------|
| R615601 | 1.5 mL conical, sterile, clear | 50/pack, 400/bag, 4000/box |
| R601001 | 15 mL, transparent rack PS     | 50/pack, 500/box           |
| R602051 | 50 mL, transparent rack PS     | 50/pack, 500/box           |

### Storage

- Store at a temperature between 15°C and 30°C.
- Keep in a dry place with relative humidity below 70%.
- Store in its original sealed packaging until ready for use.
- Do not expose to harsh chemicals or extreme temperatures that could compromise the material (polypropylene or polystyrene).